

# High Quality Surgical Scrub Sinks

Streamline™ Scrub Stations



## Built for Your OR

Streamline™ Scrub Stations are focused on your needs:

- Durable scrub sinks thermoformed of high-quality KYDEX®\* material with built-in Microban®\*\* antimicrobial protection
- Seamless large covered basin is easy to clean and reduces potential contaminants

## Hands Down - The Most Convenient Sinks

- Hands-free soap dispenser delivers 2 cc. of soap with each push of the optional foot pump accessory
- Hands-free water control – knee-actuated panel turns water on and off



Streamline single-bay scrub station



Streamline 2-bay scrub station with standard rose spray heads and optional eye wash station

## Options to Meet Surgical Staff Needs

Two different faucet head options to meet individual preferences:

- Standard “rose” spray head produces a shower-head-like pattern without aerating the water
- Optional laminar flow faucet head delivers a consistent stream of water that has a non-splash characteristic

## Reassuring Safety

- Unique eye-wash station – this faucet attachment provides convenient OR access to an emergency eye wash station:
  - when water is turned on with knee panel, one pull of the center knob releases both orange safety caps to provide instant eye flushing\*\*\*
  - when knob is not activated, faucet functions as a scrub sink
- Optimum water temperature
  - built-in thermostatically controlled mixing valve provides water temperatures of up to 115° F (46° C)



Standard rose spray head



Standard rose spray head with optional eye wash

## Scrub Sink Specifications

	1-station	2-station
Dimensions (w x d x h)	32 x 27 9/32 x 51 1/2" (813 x 692 x 1308 mm)	64 x 27 9/32 x 51 1/2" (1626 x 692 x 1308 mm)
Laminar flow faucet head	optional	optional; both faucets would have same style head
Eye wash station attachment	optional	optional; can be attached to either left or right faucet
Foot pump soap dispenser	optional	optional; recommend one for each station



Optional laminar flow head with optional eye wash

\* KYDEX® is a registered trademark of KYDEX, LLC.

\*\* Microban® is a registered trademark of Microban International.

\*\*\* Faucet-mounted eye wash stations (regardless of manufacturer) require two motions to operate (turn water on, pull knob) and thus do not meet the provisions of ANSI Z358.1-2004; these are intended as supplemental installations to other dedicated eye wash stations.

### Mixed Sources

Product group from well-managed forests, controlled sources and recycled wood or fiber



Printed with Soy based inks and water based coatings on FSC/Green Seal-certified paper.

Document # M3397EN.2010-03, Rev. A  
Printed 03/2010, 2500

©2010 STERIS Corporation.  
All rights reserved. Printed in USA.



STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 · USA  
440-354-2600 · 800-548-4873  
www.steris.com