

Harmony iQ[®] Integration Systems



Simplify Today. Embrace Tomorrow.

Simplify today.

Today's operating room is a complex environment. While patients have benefited from myriad new, life-saving technologies, OR staff are continually challenged to navigate, connect, and operate these devices. With the line of Harmony iQ Integration Systems, STERIS deployed its extensive experience in the surgical environment to create solutions that are easy to learn and use, simple to connect, and help minimize congestion.



Simplify Connections

Remove the guesswork from connecting equipment with universal ports. Cables and connectors are sturdy and durable and do not contain easily bent pins.

Simplify Control

Route the right image to the right place with an easy-to-use, intuitive touchpanel interface.




Simplify Design

Improve workflow, increase patient turn, and minimize congestion with compact systems.



Embrace tomorrow.

OR integration brings together medical imaging with audio and visual routing, two of the world's most rapidly-evolving technologies. Harmony iQ Integration Systems incorporate the latest advancements in high definition image management, bandwidth, and signal routing. Moreover, Harmony iQ Integration Systems are designed with an eye to the future, providing flexibility and scalability to ensure that your capabilities can grow as your needs evolve.



Embrace New Signals

Confidently display images from a wide variety of input sources, including all major endoscopic cameras. The end-to-end digital fiber backbone provides unrivaled image quality with flexibility for new signals in the future.

Embrace More Bandwidth

Be ready for ever-increasing bandwidth demands of evolving video technology. Get the OR integration system that is right for you today, whether your needs are big or small, complex or simple.



Embrace Advanced Capabilities

Discover upgradeable technology and virtually unlimited output destinations including video conferencing that will allow you to protect and enhance your investment.

Creating Highly Efficient Healthcare Environments

Harmony iQ® Integration Systems are available with the following options:

Surgical Grade Monitors. Outstanding resolution, color rendition, and non-glare finish. STERIS surgical monitors are light-weight and because they contain no external vents, easily cleaned.

Large-Format Wall Monitors. Exceptional HD quality. Available in 42", 47", and 55" sizes.

STERIS® Free5™ Camera System. View and record images from the surgical site from virtually any orientation and angle.

Video Conferencing. Enable one-way and two-way sharing of audio and uncompromised, HD images from the OR.

VuCapture™ Digital Recording System. Simple, integrated solution for capturing and printing surgical images and video up to 1080p.

Harmony iQ® Integration Systems

Model	2400	2800	3600
Inputs	4	8	24
Outputs	4	8	16
Weight	13.5 lbs./6.1 kg	13.5 lbs./6.1 kg	60.5 lbs/27.4 kg
Width	14.5"/36.8 cm	14.5"/ 36.8 cm	18.1"/46.0 cm
Depth	14.9"/37.8 cm	14.9"/37.8 cm	In-wall mounted: 0.8"/2.0 cm Surface mounted: 4.6"/ 11.7 cm
Height	5.8"/15.0 cm	5.8"/15.0 cm	44.7"/113.5 cm

With over 75 years in the surgical environment and the largest service organization in the industry, count on STERIS service professionals to be there when you need them - throughout the life of your equipment.



STERIS

Innovation by



Contact Your Sales Representative Today to Learn More
By Calling **1-800-548-4873**